



FOR IMMEDIATE RELEASE

For more information contact:
Cyrus Samari
651.406.9900 ext. 201

Sorna Adds Image Storage to eXpedo DICOM CD/DVD Publishing Systems

Eagan, MN, April 3, 2007 -- SToR, a new software option adding add on-line storage functionality to its highly automated intelligent eXpedo™ medical CD/DVD publishing system, was introduced by Sorna Corporation at the 2007 Convention of the American Institute of Ultrasound in Medicine held March 15-18 in New York City.

eXpedo systems are operated by intelligent, customizable DISC publishing and management software that streamlines and automates the process of creating CD/DVDs designed to contain diagnostic images (radiography, ultrasound, bone density scans, endoscopy) patient records and reports, and supplemental information.

Sorna's all-in-one XR1 and XR2 publishing systems incorporate a computer system, CD/DVD burner, a high resolution color label printer, and an LCD display with keypad controls in a compact 17" x 17" footprint. Patent-pending DISCflow™ profiles automate the production of a variety of preferred DISC formats, managing customized requests, multiple copies, and the simultaneous production of diagnostic quality and patient reference CDs.

SToR very economically adds the functionality of on-line digital image viewing from an ordinary PC or laptop using Sorna's FDA-approved DICOM viewer. This capability facilitates on-site review and discussion of radiologic or endoscopic images with patients or colleagues.

SToR also makes it possible to combine select images obtained from an external DICOM CD with new images and burn both together on a new DISC. This capability consolidates patient information contained on multiple DISCs and facilitates comparison of images using the DICOM viewer also embedded in each eXpedo DISC. SToR adds intelligence and versatility to the back-up process. Back-up can be performed on specified days and times, or scheduled when acquired data will fill a DISC.

eXpedo systems are sold through Sorna's national distribution network of established radiology dealers and distributors and OEM PACS partners.

Sorna About Sorna Corporation

Sorna is a privately held, minority owned corporation founded in October 1999. Sorna pioneered the use of 100% DICOM 3.0 compatible, automated CD/DVD equipment to provide an alternative to the use of film production for distribution and storage of radiology examinations. In addition to eXpedo CD/DVD publishing systems, its products include DICOM viewing software, DICOM DISC importing software, and digital gateways to support networked DICOM compatible hardcopy printers using Microsoft Windows™ operating software.