



## Reviewer™ Obstetrics Edition A DICOM CD/DVD Publishing System

Reviewer™ Obstetrics Edition is a versatile CD/DVD publishing system. It can easily create a DICOM DISC of diagnostic images for referral to medical specialists as well as CDs for your patients containing display-only prenatal ultrasound and fetal 3D images.

This publishing system is low-cost and easy to use. It is a high quality, manual CD/DVD burner incorporating a built-in printer that creates a full-color custom label which is printed directly onto the CD or DVD. The system is operated by software that is installed on your PC.

Like all Sorna systems, Reviewer™ Obstetrics Edition is fast and efficient. After you enter the necessary data, your professional CD will be ready in 2- 3 minutes.

**With MyBabyCD™ software, users can easily create many versatile types of CDs or DVDs.**

### Diagnostic CD or DVD

Contains DICOM 3.0 images from digital ultrasound or other DICOM sources such as – imaging equipment, archive locations or workstations. It comes with a very intuitive DICOM viewer for easy review by medical specialists or other healthcare providers.

### myBabyCD™ for Patients

Includes easy-to-display non-diagnostic ultrasound and 3D visualizations of fetal images. Additional information about your practice and medical educational resources can be included. A MyBabyCD presentation folder makes this an excellent marketing tool and a lifetime treasured memento.





## Key Features

The Reviewer Obstetrics Edition by Sorna comes with a powerful number of key features for customization and efficiency. Some of these valuable features are:

- Customizable user interface screen - modify content fields and create personal preferences.
- Search examinations residing on your PC by study and by series, with ability to select an individual image, a series, or an entire study for inclusion on the CD.
- Query any DICOM source (such as a server or workstation) by a variety of search options including patient ID, patient name, etc.,
- Export local studies to a local storage location, to an archive, or to a PACS.
- Progress bar monitors CD burning process and image search/retrieval.
- User selectable formats to add content with images. Software includes sample letter formats, description of prenatal ultrasound exams, an Internet Resource Guide for Pregnant Women, a list of best-selling books about pregnancy, and grief therapy resources.
- Support for multiple customizable labels for printing directly onto the CD in photographic quality color and detail – create different professional images for patients and referring physicians.
- The compact design and small footprint of the burner/printer occupies minimal space. Designed for professional use, it is more durable and maintenance-free than the average consumer CD/DVD burner.



myBaby CD - Table of Contents



DICOM CD/DVD - Table of Contents

### Technical Specifications

**Dimensions:** 4.75H x 6W x 10D in.

**Weight:** 7 lb. 8 oz.

**Printer Type:** Inkjet RadialPrint

**Required Media:** White/Silver Inkjet

**Data interface:** USB 2.0 high speed

**Voltage/Hz:** 110/220v 50/60Hz

**Power:** 50 watts

**Sorna Corporation**

2020 Silver Bell Road Suite 17  
Eagan, MN 55122  
USA

www.sorna.com  
1-866-767-6226