



**EMBARGOED - For Release:  
Sunday, November 26, 2006**

For more information contact:  
Cyrus Samari  
651.233.5461

**New Software for eXpedo™ Intelligent Workflow Software  
for Sorna's DICOM CD/DVD Publishing Systems Introduced at RSNA 2006**

**Eagan, MN - November 26, 2006** -- Sorna Corporation is debuting the newest version of eXpedo™, the intelligent workflow software operating Sorna's family of DICOM CD/DVD publishing systems, including the all-in-one compact eXpedo XR1 and eXpedo XR2. The new software upgrade is being demonstrated publicly for the first time in Booth 1501 South Hall at RSNA 2006.

eXpedo's patent pending software is based on the concept of DISCflow™ profiles, which define and automate the preferences of producing DICOM DISCs. eXpedo's DISCflow profiles can dramatically reduce the amount of staff time involved with DICOM DISC production and can actually increase the daily production capacity of CD/DVD burner hardware.

New enhancements include:

- Built-in DICOM CD/DVD Importer software that is fully IHE PDI Portable Media Importer Actor conformant. This new feature enables users to install the Importer on any personal computer running Windows XP.
- eXpedo Importer can also create a new DICOM DISC from the import data by sending the images and/or reports to any existing DISCflow™ profile on any eXpedo server in the network. Additionally, eXpedo Importer now has the ability to completely modify or change the contents of any DICOM field, a feature desirable when dealing with DICOM private tags.
- eXpedo now performs automatic scheduled backup of the system's database and configuration files.
- eXpedo SToR software option provides DICOM storage, Query and Retrieve functionality. Each SToR is systematically backed up to DVDs for safe, HIPAA compliant, long-term shelf storage.

**About Sorna Corporation**

Headquartered in Eagan, Minnesota, Sorna is a privately-held, minority-owned corporation founded in October 1999. It designs and manufactures 100% DICOM and IHE PDI Media profile conformant products designed to facilitate the distribution and reproduction of DICOM images and related medical data. Sorna's products include eXpedo DICOM CD/DVD publishing systems, FDA-approved DICOM viewers, the new Reviewer software, digital gateways to support networked DICOM compatible hardcopy printers with Microsoft Windows® operating platforms, and eXpedo SToR, a DICOM archive option.