



# Vertex® DISC Production Systems Comparison



**Vertex SE3**



**Vertex 15ei**



**Vertex 25ei**



**Vertex 25t**



**Vertex 22t-6**

	Vertex SE3	Vertex 15ei	Vertex 25ei	Vertex 25t	Vertex 22t-6
<b>Production Use:</b>	Light	Light/Medium	Medium	Medium/Heavy	Heavy
<b>Input Capacity:</b>	1 Bin: 20 DISCs	2 Bins: 50 DISCs ea.	2 Bins: 50 DISCs ea.	3 Bins: 50 DISCs ea.	3/4 Bins: 100 DISCs ea.
<b>Output Capacity:</b>	5 DISCs Open Bin/ 20 DISCs Closed Bin	1 Bin: 20 DISCs	10 DISCs Drawer/ 50 DISCs Bin	50 DISCs Internal Bin/ 50 DISCs External Bin	100 DISCs Internal Bin/ 100 DISCs External Bin
<b>Recorder(s):</b>	1 Recorder: Burns CD-R or DVD-R	1 Recorder: Burns CD-R & DVD-R	2 Recorders: Burns CD-R & DVD-R	2 Recorders: Burns CD-R & DVD-R	2 Recorders: Burns CD-R & DVD-R
<b>Printing Technology:</b>	High-Resolution Inkjet	High-Resolution Inkjet	High-Resolution Inkjet	High-Resolution Black or Color Thermal	High-Resolution Black or Color Thermal
<b>Media:</b>	Standard Inkjet CD-R or DVD-R	Standard Inkjet CD-R & DVD-R	Standard Inkjet CD-R & DVD-R	Thermal Transfer Certified CD-R & DVD-R	Thermal Transfer Certified CD-R & DVD-R
<b>Dimensions:</b>	11.4" x 3.7" x 12.3"	11.4" x 13.7" x 6.8"	11.4" x 13.7" x 12.3"	22.8" x 15.1" x 21.5"	14.2" x 13.7" x 6.8"
<b>Highlights:</b>	<b>Low-cost</b> , efficient, & reliable DICOM DISC publishing system with custom DISC labeling.	<b>Affordable</b> , high- capacity, & efficient at a compact size.	High-performance, high- capacity on a <b>moderate budget</b> with a pullout output bin.	High-volume system for DISC production. <b>Best price</b> for performance thermal printing.	Provides <b>maximum performance</b> , reliability, capacity, & thermal print speed.