

Budget friendly and easy to use software for exporting studies to local and removable drives. Add a viewer and choose to encrypt the data protecting patient PHI.

With its ability to copy to removable media as well as local and network drives, Export is best used for small imaging centers, or as a backup to an automated CD burner. It can also be utilized to anonymize data for teaching purposes. Add Vertex SHARE to enable cloud brokered image sharing.

Create standards compliant DISCs: Utilizing Vertex's intuitive interface, users are able to burn studies to local CD-R, DVD-R, and BD-R drives – fully conformant to DICOM and IHE standards.

Export studies to removable Media: DICOM studies can be copied to USB and other removable media for distribution to patients, healthcare facilities, and referring physicians.

PHI Security: Choose only the studies and series you wish to import. Eliminate the need to import non-essential studies – saving time, storage space, and reducing liability.

Copy studies to drives: Share studies directly with authorized professionals within your facility using local and network drives.

Embedded Viewer: Developed in-house, our intelligent and easy-to-use embedded viewing software improves your referring physician and patient experience.

PHI Security: Encrypt the data to comply with DICOM and FIPS 140-2 requirements. The required password can be generated by four different methods, depending on preferred workflow (on-demand, pre-set, DICOM tag, random).


One-Click De-Identification: Quickly and easily create teaching files and/or research projects

Track all activity: Data exported from Vertex are tracked within its database for HIPAA compliance

Unlimited client License: Any number of authorized users, groups, or workstations can utilize Vertex Export without running into licensing limitations or extra fees

Client Access Control: For security and regulatory purposes, individual users, groups, or workstations can be configured to have access to incoming devices and export destinations

Workflow simplifications: Minimize operator responsibility and reduce human errors with predefined routing rules.



Vertex Platform

Medical imaging made easy- the most affordable, scalable medical data systems



Cloud-brokered medical data sharing



Automated DICOM CD burners plus CD importing



Post-process and convert exiting video files to DICOM



Compact, economical video recorder integrated with Vertex Video



Add Paper Film, Digital Images and Documents to PACS- scan paper or film and include images and documents in studies



Reconcile Incoming DICOM and Proprietary Studies- Import DICOM and proprietary data from CDs and removable media



Manually send to local CD burner, removable media, and drives



Securely send studies using DICOM TLS with no VPN



Scheduled, Intelligent routing to any device

REVIEWER

Manage internally and externally generated DICOM and non-DICOM medical images within your role and workflow with Reviewer.

Reviewer CD

View and browse files on every CD (DICOM and proprietary) with this easy-to-use, FREE DICOM viewer software

Reviewer MD

Burn, view, store, import, and share studies from a variety of sources; send to a printer or create an encrypted DISC on the fly

Reviewer Admin

View, import, burn, share, scan, re-identify, encrypt, and MORE; do it all with this Master Edition



info@sorna.com • 651.406.9900