

vertex[®]

Version 3.8

IHE Integration Statement



Contents

Introduction	3
Revision History	3
Legal Notices	3
Copyright Notice	3
Trademark Acknowledgements	3
Contact Address	3

Introduction

Revision History

Revision	Author	Comments	Date
1	Dion Barrier	Created	04/20/2016
2	Dion Barrier	Reviewed, updated	04/29/2020

Legal Notices

Copyright Notice

©2020 Sorna Corporation. All rights reserved.

- No part of this document may be copied or reprinted, in whole or part, without written permission.
- The contents of this document are subject to change without notice or legal obligation.

Trademark Acknowledgements

Vertex is a registered trademark of Sorna Corporation. This document may include trademarks or registered trademarks of other companies.

Contact Address

Sorna Corporation
2020 Silver Bell Road
Suite 17
Eagan, MN 55122

<i>IHE Integration Statement</i>			
Vendor	Product Name	Version	Date
Sorna Corporation	Vertex	v3.8	04/29/2020
This product is indented to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below.			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Import Reconciliation Workflow (IRWF.b)	Importer	Scheduled Import - Uses MWL	
Portable Data for Imaging (PDI)	Portable Media Creator		
	Portable Media Importer		
	Report Reader		
Internet address for vendor's IHE information: www.sorna.com			
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	