

vertex®

Enterprise Imaging - Simplified

Software Release

Version: 3.9

Release Date: 10-01-2021

New Features:

Introducing Enterprise CD. Expanding the standard Vertex CD set up, Enterprise CD allows you to:

- Securely and conveniently setup, configure, and access multiple CD burners.
- Achieve enterprise-wide CD burning by physically separating CD burners from the server (run Vertex on a VM)
- Single database tracking all activity across the enterprise
- User login provides access to all permitted autoloaders
- Run multiple Vertex applications from a central Vertex platform installation (e.g., SHARE, Scan, Route, ...)

Spotlight Enhancements:

- Added support for the cost effective 24t autoloader
- Added support for Maestro (printable USB autoloader)
- Improved TWAIN scanner management
- Improved Windows / Active Directory Authentication Options

Contact your Sorna Representative to learn more about the new features and additional enhancements now available. See if your site qualifies for the free upgrade; software upgrades are free to any site under an Advanced Exchange Maintenance contract.