







## Robust CD/DVD/USB burning solutions to meet the needs of your healthcare facility.

**Packages include:** Vertex CD+ software running on customized PC, burner, first year Advanced Exchange Maintenance, and starter supply kit. Enterprise configuration available.

Vertex SE3	Highlights	Production	Input Capacity	Output Capacity	Recorder(s) and Media:	Printing Technology:	Dimensions (Approx. H x W x D)
	Low-cost & efficient DICOM DISC publishing system with custom labeling	Light	1 Bin: 20 DISCs	Open Bin: 5 DISCs Closed Bin: 20 DISCs	1 Recorder: Burns CD-R or DVD-R	High-Res Inkjet on Inkjet CD-R & DVD-R	7" x 15" x 15"
<b>Vertex 15ei</b>							
	Affordable, high capacity, & efficient in a compact size	Light/ Medium	2 Bins: 50 DISCs ea.	1 Bin: 20 DISCs	1 Recorder: Burns CD-R & DVD-R	High-Res Inkjet on Inkjet CD-R & DVD-R	13" x 15" x 18"
<b>Vertex 25ei</b>							
	High-performance, high-capacity on a moderate budget - pullout output bin	Medium	2 Bins: 50 DISCs ea.	Drawer: 10 DISCs Bin: 50 DISCs	2 Recorders: Burns CD-R & DVD-R	High-Res Inkjet on Inkjet CD-R & DVD-R	15" x 15" x 20"
<b>Vertex 24t</b>							
	Compact and high-performance, best priced thermal DICOM CD burner	Medium	3 Bins: 50 DISCs ea.	External Bin: 5 DISCs <i>(optional)</i> External Bin: 50 DISCs	2 Recorders: Burns CD-R & DVD-R	High-Res. Black or Color Thermal Transfer CD-R & DVD-R	19" x 15" x 20"
<b>Vertex 25te</b>							
	High-volume system for DISC production with performance thermal printing	Medium/ Heavy	3 Bins: 50 DISCs ea.	50 DISCs Internal Bin/ 50 DISCs External Bin	2 Recorders: Burns CD-R & DVD-R	High-Res. Black or Color Thermal Transfer CD-R & DVD-R	26" x 16" x 23"
<b>Vertex RX4</b>							
	Compact system for quickly recording large amounts of data on USB flash drive	Light/ Manual	4 USB Slots	NA	4 USB Slots: Multi-use or Single-use USB sticks (16/64GB)	Paper Label / Black	7" x 18" x 10"