

vertex[®]

Enterprise Imaging - Simplified

Software Release

Version: 3.10

Release Date: 06-13-2023

New Vertex Platform Products:

Vertex CD is now Vertex CD+

- Automated DICOM CD burning, USB production; and importing.

Vertex Upload

- Vertex Upload is a cloud-brokered solution for patients and providers to securely send medical data.
- Email a link and a separate Access Code.
- No account set up and no cloud storage.
- View received DICOM and non-DICOM datasets in your on-site Vertex Inbox.

Vertex RX4 USB Production System

- Network-attached, four-slot USB production system with label printer.
- Integrated with Vertex CD+.
- Fast and secure distribution of large datasets on USB media.

Spotlight Enhancements:

- Philips Intellispace PACS DICOM Report Interface.
- Audit records for administrative and PHI query activities.
- Series level DICOM query/retrieve (if supported by the archive).
- DICOM filter/routing rule now supports string "contains."
- Override DICOMization template study date/time.
- Security enhancements.

Contact your Sorna Representative to learn more about the new features and additional enhancements now available. See if your site qualifies for the free upgrade; software upgrades are free to any site under an Advanced Exchange Maintenance contract.