

The vertex logo consists of the word "vertex" in a bold, black, sans-serif font. The letter "v" is stylized with a green triangle pointing downwards from its top left corner and a white horizontal bar across its middle.

Version 3.10

IHE Integration Statement

Contents

| | |
|--|---|
| Introduction | 3 |
| Revision History..... | 3 |
| Legal Notices | 3 |
| Copyright Notice..... | 3 |
| Trademark Acknowledgements..... | 3 |
| Contact Address | 3 |
| <i>IHE Integration Statement</i> | 4 |

Introduction

Revision History

| Revision | Author | Comments | Date |
|----------|--------------|-------------------|------------|
| 1 | Dion Barrier | Created | 04/20/2016 |
| 2 | Dion Barrier | Reviewed, updated | 04/29/2020 |
| 3 | Dion Barrier | Reviewed, updated | 11/01/2023 |

Legal Notices

Copyright Notice

©2023 Sorna Corporation. All rights reserved.

- No part of this document may be copied or reprinted, in whole or part, without written permission.
- The contents of this document are subject to change without notice or legal obligation.

Trademark Acknowledgements

Vertex is a registered trademark of Sorna Corporation. This document may include trademarks or registered trademarks of other companies.

Contact Address

Sorna Corporation
4101 Nicols Road
Suite 100
Eagan, MN 55122

| <i>IHE Integration Statement</i> | | | |
|---|--|---|-----------------------------|
| Vendor | Product Name | Version | Date |
| Sorna Corporation | Vertex | v3.10 | 11/01/2023 |
| This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below. | | | |
| Integration Profiles Implemented | Actors Implemented | | Options Implemented |
| Import Reconciliation Workflow (IRWF.b) | Importer | | Scheduled Import – Uses MWL |
| Portable Data for Imaging (PDI) | Portable Media Creator | | |
| | Portable Media Importer | | |
| | Report Reader | | |
| Internet address for vendor’s IHE information: www.sorna.com | | | |
| Links to general information on IHE | | | |
| In North America: www.rsna.org/IHE | In Europe: www.ihe-europe.org | In Japan: www.jira-net.or.jp/ihe-j | |