

**Tired of installing scanning software on every workstation? Using multiple applications to convert digital documents? Vertex Scan is just a click away every time.**

Customize the workflow of adding paper, film, digital images and documents to existing studies or creating new ones. Affordability, flexibility, and ease-of-use are the hallmarks of Vertex Scan.

**Scan All Papers:** Add consent forms, modality worksheets, insurance cards, IDs, and other documents to studies.

**Create Studies from Film:** Digitize films from prior studies to store in PACS or share with healthcare facilities and patients.

**Manage Visible Light Images:** DICOMize popular image formats (JPEG, BMP, PNG, TIFF, GIF, ...) for inclusion in PACS and VNA.

**Add Digital Documents:** Office documents such as PDF, Word, Excel, Text, HTML, RTF, etc. can be DICOMized and seamlessly added into studies.

**Multiple Modes of Operation:** Initiated from PACS workstation, papers, film, documents, and images can automatically be added to studies. Users can also drag-and-drop any number of objects into an existing study or create a new one.

**Flexible DICOMization Options:** Select pagination format for document conversion. PDFs can be stored as embedded or secondary capture objects.

**Create Unique Studies:** Use Vertex to manage automatic creation of unique accession numbers and/or connect to Modality Worklist Server.

**Consistent DICOMization:** Configure department information, and create pre-defined study and series descriptions, ensuring data integrity when creating new studies.

**HIS Friendly:** Add procedure codes to pre-defined study descriptions to ensure proper billing and reimbursement

**Client Access Control:** For security and regulatory purposes, individual users, groups, or workstations can be configured to have access to incoming devices and destinations.

**License All Scanners:** Select the number of scanners to license.

