

Vertex® 3.11.4.0 Marketing Release Notes

Major enhancements and feature additions implemented for Vertex 3.11.4.0 are:

Event Notification Templates

- Customizable email and SMS event notification templates
- Merge text fields
- Support for Vertex SHARE and Upload

Automatic DICOM Import Reconciliation

- Automates the selection and querying of the reconciliation device based on prior configuration
- Pre-fills search queries using relevant context from the source data
- Streamlines the normalization process against in-facility sources such as PACS or a Modality Worklist by reducing manual input
- Supports all Vertex source devices (removable media, local and network drive(s), cloud, ...)

Vertex Workstation Enhancements

- Searchable address book for cloud providers
- Trial/term limited license status User Interface
- Patient DOB added to worklist query form
- Harmonized device query fields in the GUI

General Enhancements

- Support for RX-4 international character sets
- Flat-folder client file system export option
- Media template determines transfer-syntax export behavior
 - support automatic compression of data export to autoloaders
- 25ei maintenance box supply low/out event notification
- License expiration warning/notice event notification
- Dataset size-based unattended profile routing rules